



Isolator Mounting Hole Dimensions

WIDTH (W)	LENGTH (L)	DIMENSION (A)
in. (mm)	ft. (mm)	in. (mm)
36.0 (900)	20.0 (6000)	53.0 (1346)
48.0 (1200)	18.0 (5400)	47.0 (1194)
59.1 (1500)	16.0 (4800)	42.0 (1067)
	14.0 (4200)	37.0 (940)
	12.0 (3600)	32.0 (813)
	10.0 (3000)	27.0 (686)
	8.0 (2400)	22.0 (559)
	6.0 (1800)	16.0 (406)